

# The effect of intermolecular forces and entropy on viscosity, vapor pressure and evaporation in hydrocarbon oils

Tom Karis, Raj Shah, Andrew Zheng



## Motivation

Hydrocarbon base oil provide the starting point for lubricant formulations. The fundamental problem addressed in this work is the relationship between viscosity and vapor pressure of base oil.


Needs for low viscosity/low vapor pressure oils.

- Increase efficiency of motors, engines, and drivetrains
- Lower friction and evaporation loss

## Problem:

- Low viscosity -> leads to high vapor pressures


## Vapor Pressure vs. Viscosity for Two Base Oils




Acronym	Molecular Weight (g/mol)	Structure
PAO	240	<chem>CCCCCCCCCCCCCCCC</chem>
DOS	427	<chem>CCCCCCCCCCCCCCCCC(=O)OCCCCCCCCCCCCCCCC</chem>

## “Arrhenius” plots for Vapor Pressure and Viscosity

The slope and intercept determine the vapor pressure and viscosity




Acronym	Molecular Weight (g/mol)	Structure
DOA	371	<chem>CCCCCCCCC(=O)OCCCCCCCCC(=O)OCCCCCCCCC</chem>
DIA	427	<chem>CCCCCCCCC(=O)OCCCCCCCCC(=O)OCCCCCCCCC(C)(C)C</chem>



## References


1. T. E. Karis & H. S. Nagaraj (2000) Evaporation and Flow Properties of Several Hydrocarbon Oils, TRIBOLOGY TRANSACTIONS, 43:4, 758-766, DOI: 10.1080/10402000008982405
2. R.H.Ewell & H.Eyring (1937) Theory of Viscosity of Liquids as a Function of Temperature and Pressure, Journal of Chemical Physics, Princeton, NJ
3. S.Q.A. Rizvi, Petroleum-Derived Hydrocarbon Base Oils, Fuels and Lubrications Handbook ©2016, pg 292
4. J.H. Hildebrande(1939), Liquid Structure and Entropy of Vaporization, Journal of Chemical Physics, Berkeley, CA
5. W. Kauzmann & H. Eyring (1940) The Viscous Flow of Large Molecules, J.AM Chemistry Society Vol. 62. Princeton, NJ

## Evaporation: Forms a Hole




vaporization entropy  $\Delta S_{vap} \approx S_{vap} - S_{liq}$

## Flow: Changes Position



flow activation entropy  $\Delta S_{vis} = \Delta S_{vis}^{trans} + \Delta S_{vis}^{rot}$

## Flow Activation Rotational Entropy $\Delta S_{vis}^{rot}$ vs. Shape Factor $n$




Molecules with a higher vaporization energy have a higher viscosity because they gain less or lose more rotational entropy as they transition to a hole in the liquid.

The slope and intercept of the Arrhenius plot are fixed by the activation entropy  $\Delta S$  and energy  $\Delta H = \Delta E + \Delta(PV)$  for vaporization and for flow.

$P_{vap}$  vapor pressure,  $P_{atm}$  atmospheric pressure,  $\mu$  kinematic viscosity,  $M$  molecular weight,  $N$  Avogadro's number,  $h_p$

Combining the vapor pressure and viscosity activation energies, the molecular shape factor is provided by the vaporization  $\Delta E_{vap}$  and flow  $\Delta E_{vis}$  activation energies (determined by the number and type of atoms and oil polarity)



## Conclusion

Lubricant base oils with lower viscosity and vapor pressure should be designed to gain more entropy as they translate between the neighboring molecules in the liquid. This can most likely be accomplished by designing a molecular structure with a greater flow rotational activation entropy and lower vaporization entropy.